

DJ Topholm

DNK 302046

Kappa casein: BB
DOB.: 21/11/2000
Proven

aAa: 615243
Toptot-HH x JAS Byg x SKÆ Krig



IMPORT



NTM -11

RECOMMENDED FOR

- HOOF HEALTH

Dam average (305 days) Data not available
Est. dtr. average (305 days) 7310 kg 6.29 % F / 452 kgF 4.31 % P / 315 kgP

NTM July 2024

PRODUCTION & EFFICIENCY NO. OF DTRS 3510 REL 99%

		90	100	110	120
Production	91	[Bar chart showing 91 vs 90-120]			
Milk kg	88	[Bar chart showing 88 vs 90-120]			
Protein kg	83	[Bar chart showing 83 vs 90-120]			
Protein %	100	[Bar chart showing 100 vs 90-120]			
Fat, kg	97	[Bar chart showing 97 vs 90-120]			
Fat, %	112	[Bar chart showing 112 vs 90-120]			
Growth	109	[Bar chart showing 109 vs 90-120]			
Saved feed	104	[Bar chart showing 104 vs 90-120]			
Persistency	94	[Bar chart showing 94 vs 90-120]			

FUNCTIONAL NO. OF DTRS 2868 REL 99%

		90	100	110	120
Daughter fertility	88	[Bar chart showing 88 vs 90-120]			
Calving sire	94	[Bar chart showing 94 vs 90-120]			
Calving maternal	84	[Bar chart showing 84 vs 90-120]			
Udder health	108	[Bar chart showing 108 vs 90-120]			
General health	94	[Bar chart showing 94 vs 90-120]			
Longevity	97	[Bar chart showing 97 vs 90-120]			
Hoof health	111	[Bar chart showing 111 vs 90-120]			
Youngstock survival	99	[Bar chart showing 99 vs 90-120]			
Milkability	87	[Bar chart showing 87 vs 90-120]			
Temperament	104	[Bar chart showing 104 vs 90-120]			

OVERALL CONFORMATION NO. OF DTRS 2118 REL 99%

		90	100	110	120
Frame	80	[Bar chart showing 80 vs 90-120]			
Feet & legs	104	[Bar chart showing 104 vs 90-120]			
Udder	85	[Bar chart showing 85 vs 90-120]			

CONFORMATION

		80	90	100	110	120	
Stature	Short	[Bar chart showing Short vs 80-120]				Tall	88
Rib structure	Coarse	[Bar chart showing Coarse vs 80-120]				Angular	87
Chest width	Narrow	[Bar chart showing Narrow vs 80-120]				Wide	94
Body depth	Shallow	[Bar chart showing Shallow vs 80-120]				Deep	94
Rump width	Narrow pins	[Bar chart showing Narrow pins vs 80-120]				Wide pins	76
Rump angle	High pins	[Bar chart showing High pins vs 80-120]				Low pins	103
Top line	Weak	[Bar chart showing Weak vs 80-120]				Strong	91
Rear legs, side view	Straight	[Bar chart showing Straight vs 80-120]				Sickled	101
Rear legs, rear view	Toes out	[Bar chart showing Toes out vs 80-120]				Parallel	97
Foot angle	Low	[Bar chart showing Low vs 80-120]				Steep	116
Bone quality	Coarse	[Bar chart showing Coarse vs 80-120]				Fine and thin	101
Hock quality	Filled	[Bar chart showing Filled vs 80-120]				Dry	95
Fore udder attachme...	Loose	[Bar chart showing Loose vs 80-120]				Strong	85
Rear udder height	Low	[Bar chart showing Low vs 80-120]				High	88
Rear udder width	Narrow	[Bar chart showing Narrow vs 80-120]				Wide	87
Udder support	Weak	[Bar chart showing Weak vs 80-120]				Strong	91
Udder depth	Deep	[Bar chart showing Deep vs 80-120]				High	84
Udder balance	Deep rear	[Bar chart showing Deep rear vs 80-120]				Deep front	93
Teat length	Short	[Bar chart showing Short vs 80-120]				Long	91
Teat thickness	Thin	[Bar chart showing Thin vs 80-120]				Thick	99
Teat placement (front)	Wide	[Bar chart showing Wide vs 80-120]				Close	78
Teat placement (back)	Wide	[Bar chart showing Wide vs 80-120]				Close	78